


RADON

MEDICAL IMAGING



Sorna Vertex Share

radonmedicalimaging.com

We focus on people, innovative products,
and service excellence

WV Office 866-723-6698

VA Office 800-722-1991

What is it?


Vertex SHARE is a cloud-brokered solution for securely sharing medical data with providers and patients.

Share

Use Vertex to collate and share DICOM studies, paper, film, office documents, visible light images, video, and other digital content. Selected recipients will receive an email notification with an Access Link to review and download the shared data. No account setup required, just enter the Access Code provided when the SHARE was created.

Download

No special software required, just an ordinary web browser. Recipients are provided with a choice of download formats—either a ZIP archive, or a virtual CD/DVD (ISO image) file—as well as the option to include viewing software. Large datasets such as those associated with Release of Information requests are handled with ease thanks to multi-download spanning feature of SHARE.



Secure

End-to-end encryption

Data in transit is always encrypted. Network transport utilizes TLS v1.2. Static data encrypted with AES-CBC 256bit.

No in-bound firewall exceptions

Vertex SHARE utilizes IT-friendly outbound connections only.

Least-privilege access

Actors at each step of the process between share notification through to download are restricted to only those resources required for their part.

In Control

As the sharing facility, you are in control of your data. You control:

What is shared ... You choose what data to share.

Where its cached ... All medical data is retained locally – not the cloud.

Who receives it ... Vertex SHARE supports both established relationships as well as ad-hoc requests.

How long ... You determine your own retention policy.

In the Know

Vertex maintains a complete, searchable history of SHARE activity. Know when, what and with whom it was shared; and just as important, when and if it were downloaded.

vertex[®]SHARE

Securely Share Medical Data with Patients and Medical Professionals



How do I share?

1. Using Vertex; combine, edit, and package the data to be shared.
2. Select Vertex Share, enter email address(es) or select from the list.
3. Choose any additional destinations; click submit.

How do they access?

1. Receive the email with a time-sensitive link to the data package.
2. Click on the link and input password to gain access to the content.
3. Select data format and viewer; click Download.

- » Medical data is retained locally and streamed on demand through the Vertex Share cloud.
- » Streaming minimizes upload bandwidth, eliminates cloud storage costs and significantly reduces security risks.
- » DICOM studies, paper, and other digital content are combined by Vertex into a single package ready for sharing.
- » Packaged data can be sent to multiple Destinations (CD, USB, ...) under the same job – eliminating duplication of effort.
- » No account set up is necessary for recipients – simplify their process, save time and money.

Contact Sorna today for a customized demonstration of Vertex Share.