


RADON

MEDICAL IMAGING



Sorna 15ei CD Burner

radonmedicalimaging.com

We focus on people, innovative products,
and service excellence

WV Office 866-723-6698

VA Office 800-722-1991



Harmony. Innovation. Simplicity.



Vertex 15ei

High-quality output: *Micro Piezo printhead*

Economical printing: *Individual inks system*

Reliable operation: *Dust-proof body*

Accurate grip: *Patented robotic arm*

Easy to use: *Front operation*

Production Use: Light / Medium

Input Capacity: 2 Bins: 50 DISCs ea.

Output Capacity: 1 Bin: 20 DISCs

Recorder(s): 1 Recorders: Burns CD-R or DVD-R

Printing Technology: High-Resolution Inkjet

Media: Standard Inkjet CD-R & DVD-R

Dimensions: 12.76" x 14.84" x 18.31"

Highlights: Affordable, high-capacity, and efficient at a compact size



Visit www.sorna.com

Patent: www.sorna.com/patents

*Image not to scale.

Enterprise Imaging - Simplified.

Sorna Corporation's industry leading Vertex[®] Platform software is a single application that can combine video, images and data from across the healthcare enterprise, convert them to DICOM for storage in PACS/VNA and share them via CD/DVD/cloud/USB and more.

vertex[®]
powered CD / DVD
burning solutions



Vertex SE3



Vertex 15ei



Vertex 25ei



Vertex 24t



Vertex 25te

Systems:	Vertex SE3	Vertex 15ei	Vertex 25ei	Vertex 24t	Vertex 25te
Highlights:	Low-cost, efficient, & reliable DICOM DISC publishing system with custom DISC labeling	Affordable, high-capacity, & efficient in a compact size	High-performance, high-capacity on a moderate budget with a pullout output bin	Compact and high-performance, best priced thermal printing DICOM CD burner	Reliable, high-volume system for DISC production with performance thermal printing
Production:	Light	Light / Medium	Medium	Medium	Medium / Heavy
Input Capacity:	1 Bin: 20 DISCs	2 Bins: 50 DISCs ea.	2 Bins: 50 DISCs ea.	3 Bins: 50 DISCs ea.	3 Bins: 50 DISCs ea.
Output Capacity:	5 DISCs Open Bin/ 20 DISCs Closed Bin	1 Bin: 20 DISCs	10 DISCs Drawer/ 50 DISCs Bin	5 DISCs External Bin (optional 50 DISCs External Bin)	50 DISCs Internal Bin/ 50 DISCs External Bin
Recorder(s) and Media:	1 Recorder: Burns CD-R <u>or</u> DVD-R	1 Recorder: Burns CD-R & DVD-R	2 Recorders: Burns CD-R & DVD-R	2 Recorders: Burns CD-R & DVD-R	2 Recorders: Burns CD-R & DVD-R
Printing Technology:	High-Resolution Inkjet Standard Inkjet CD-R & DVD-R	High-Resolution Inkjet Standard Inkjet CD-R & DVD-R	High-Resolution Inkjet Standard Inkjet CD-R & DVD-R	High-Resolution Black or Color Thermal Transfer CD-R & DVD-R	High-Resolution Black or Color Thermal Transfer CD-R & DVD-R
Dimensions (H x W x D):	7" x 15" x 14.75"	12.7" x 14.8" x 18.3"	13.7" x 14.8" x 20.2"	19" x 14.8" x 19.5"	25.5" x 16" x 23.25"

Enhance your facility's efficiency with Vertex's cost-effective and versatile array of DICOM CD Burner hardware systems for any production level and budget.

Vertex CD is the latest automated DICOM CD burner software technologies

IT pioneered by Sorna to simplify enterprise imaging. The Vertex family of recording systems targets every production level and budget. Select from a large array of affordable inkjet or thermal CD / DVD burning solutions to meet the needs of each department within your healthcare facility.

Easy to Learn and Use

With a single screen and drag-and-drop functionality, Vertex significantly simplifies the user experience and accelerates operator training.

Import DICOM and Proprietary DISCs

Transparently import poorly formatted and proprietary DISCs into your PACS-eliminating the need to re-scan a patient or request another DISC.

Attended Workflow

Operators can complete multiple tasks simultaneously for maximum production.

Unattended Workflow

Customizable data profiles-automatically share studies with multiple destinations.

Installation-Free Client Deployment

Click-once deployment technology minimizes IT requirements and streamlines the workstation installation and maintenance process.

Integration with Existing Hardware

Maximize your current hardware investment through a seamless integration with all industry-leading auto loaders.

PACS Integration and Support

Developed to easily connect with your PACS in minutes, reducing installation cost.

Virtually Limitless Data Profiles

Vertex handles numerous jobs concurrently, substantially reducing production time and enhancing efficiency.

DISC-Embedded Viewing

Developed in-house, our intelligent and easy-to-use embedded viewing software improves your referring physician and patient experience.

One-Click De-Identification

Quickly and easily create teaching files and/or research projects.

Encryption

Four Methods Secure DICOM and FIPS 140-2 conformance.

Transaction Tracking

Track all transactions for HIPAA compliance.



Vertex CD Burner System Package

- 1- Vertex Server and Client software licenses per system and DICOM Import of medical images in standard and proprietary formats from CD, USB, and drives.
- User Query & Build, Automated DICOM processing, and popular PACS interfaces are also provided.
- Produced CDs include a DICOM Viewer for PCs, a DICOM viewer for Mac users, and an installable free DICOM viewer to view all medical CDs.
- First year Advanced Exchange maintenance
- A customized-for-user PC with monitor, keyboard, mouse, and all cabling
- 200 CDs/DVDs with holders and ink cartridge