

# MED-RX Contrast Injector Tubing

Radon offers significant savings on disposable tubing used primarily on MedRad injectors.

Contrast Delivery - Cross Reference										
Product Description	MED-RX®		MEDRAD		Coeur		Covidien		Merit	
	Code	Case Qty.	Code	Case Qty.	Code	Case Qty.	Code	Case Qty.	Code	Case Qty.
<b>Fluid Transfer (Fill) Sets</b>										
9" Transfer (Fill) Set, Female Luer, Vented Spike, Pinch Clamp	10-1302	50			C405-3041	100				
21" Transfer (Fill) Set, Swabbable Valve, Vented Spike	10-1304	100			C405-3200	100			TS220	50
21" Transfer (Fill) Set, Swabbable Valve, Ball Check Spike	10-1306	50	SVTS*	50						
21" Transfer (Fill) Set, Swabbable Valve (Large Wings), Ball Check Spike	10-1308	100	SVTS*	50	C405-3208	100			TS120	50
<b>Y-Tubing Sets</b>										
18" Y-Tubing Set, Dual Check Valve, 400 PSI Rated	10-1253	50								
70" Y-Tubing Set, Dual Check Valve, Coiled Line, 400 PSI Rated	10-1257	100			C405-2155	100	844012	50		
70" Y-Tubing Set, Single Check Valve, Coiled Line, 400 PSI Rated	10-1258	100			C405-2604	100	844011	50		
18" Y-Tubing Set, Single Check Valve, 400 PSI Rated	10-1259	100								
70" Y-Tubing Set, No Check Valve, Coiled Line, 400 PSI Rated	10-1260	100			C405-1549	100	844010	50		
<b>Coiled Line</b>										
60" Low Pressure Coiled Line, Male & Female Connectors, 400 PSI Rated	10-1309	100	LPDCT-160	50	C405-0060	100	601195	100		
		Rated 400 psi				Rated 300 psi		Rated 400 psi		

## Budget Pricing (sample of products):

Description	Qty.	Unit Price	Per Set pricing
10-1308	Per Case (100 sets)	\$319	Per set price \$3.19
21" Transfer (Fill) Set, Swabbable Valve (Large Wings), Ball Check Spike			
10-1306	Per Case (50 sets)	\$148	Per set price \$2.96
21" Transfer (Fill) Set, Swabbable Valve, Ball Check Spike			
10-1304	Per Case (100 sets)	\$243	Per set price \$2.43
Transfer Set with Swabbable valve			
10-1302	Per Case (50 sets)	\$69	Per set price \$1.38
Transfer Set 9" with female Luer, C clamp with vented spike			
10-1309	Per Case (100 sets)	\$155	Per set price \$1.55
60" Coiled Line, Rated to 400 PSI			

# MED-RX Contrast Management Products


- FDA Approved
- Specs meet or exceed OEM Manufacturer specs
- Easier connection (improved flexibility over OEM design)
- Saving of 25% or more off OEM pricing

## 60" Coiled Lines

"10-1309"



- Non-DEHP PVC tubing
- Male Luer is vented in the center for easy Purging of the line
  - Rated to 400PSI
- Packaged sterile, 100/case




- Y sets also available with coiled lines
- Available with and without check valves

# MED-RX Contrast Management Products (continued)


## Contrast Transfer Sets

### "10-1302"




- Transfer Set 9" with female Luer, C clamp with vented spike
  - Non-DEHP PVC tubing
  - Packaged sterile, 50/case

### "10-1308"



### "10-1304"



- Transfer Sets with swabbable valve
- 10-1308 has vented spike with ball check valve
  - Non-DEHP PVC tubing
  - Packaged sterile, 100/case



### "10-1306"

21" Transfer (Fill) Set,  
Swabbable Valve, Ball Check  
Spike (300 PSI)

- Specs meet or exceed OEM Manufacturer specs
- Easier connection (improved flexibility over OEM design)
- Saving of 25% or more off OEM pricing