

# Weight Bearing Protection



*Our new Weight Bearing Protective Cover has been designed to accommodate Digital DR Panels and CR Cassettes.*

*Our covers are strong yet lightweight, offering a high degree of protection during weight bearing procedures.*

*Our latest cover is constructed of a high density polyethylene frame and a carbon fiber top, offering protection up to 1300 lbs. Our easy lift handle makes removal and carrying a breeze.*

*The Weight Bearing Cover has been tested by an accredited independent laboratory.*



# Weight Bearing Protection



Smoked Polycarbonate



Carbon Fiber

- Available in Polycarbonate or Carbon Fiber
- Rated up to 1300 lbs.
- Easy lift handle
- Non-skid rubber feet
- Multiple sizes available
- Printed target area

*This product is perfect for bariatric, equine, and podiatry procedures.*



Easy Lift Handle

Protect your *investment* use a **COVER!**



**RC Imaging**  
your radiology resource

800.528.6311 585.392.4336 585.392.5575  
50 Old Hojack Lane, Hilton, New York 14468  
[www.rcimaging.com](http://www.rcimaging.com)

